



# MD **Gekko**<sup>TM</sup> and SiC Foil

## Easy applying and removing



*Fixation is child's play!*

**The new MD-Gekko with the unique surface makes applying and removing all types of consumables easier than ever before.**

The only requirement to the consumables applied is that the backing has to be flat and glossy, like with PET backed consumables. The grinding or polishing consumable clings safely to the MD-Gekko during preparation and after use the consumable can be peeled off, without effort leaving no trace of residue on the surface of MD-Gekko. This also means they can be reused until completely worn.

Traditionally, consumables have been fixed using glue or double adhesive tape, often making them difficult to apply, difficult to remove and impossible to reuse after having been removed. With MD-Gekko this is all history.

The unique surface allows even self-adhesive consumables to be attached and removed easily, without leaving any trace of adhesive.

## No wrinkling or wrapping

Together with MD-Gekko we also introduce a new SiC Foil, to be used on MD-Gekko. The new SiC Foil is based on a glossy PET foil, the perfect match for MD-Gekko. Water does not affect the shape of the PET foil at all. Thus it does not wrinkle or wrap as traditional SiC Paper, but stays flat and ready for immediate or later re-use.

Changing between the 11 available grit sizes, from # 80 for coarse grinding to # 4000 for fine grinding, is done in the shortest possible time.

It is not necessary to apply water to a disc prior to use or carefully apply paper on a self adhesive disc, just place the SiC Foil on MD-Gekko and you are ready to go.

**Grinding on SiC has never been easier!**



Traditionally, consumables have been fixed using glue or double adhesive tape, making them often difficult to apply, difficult to remove and impossible to reuse after having been removed. With MD-Gekko this is all history. The unique surface allows even self-adhesive consumables to be attached and removed easily, without leaving any trace of adhesive.

**MD-Gekko can be used for a wide variety of consumables again and again and again...**

# NEW

**Struers**

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk  
www.struers.com